



國立中正大學

National Chung Cheng University



積極創新 修德澤人



Prof. Joyce S. Yu

Department of Chemistry and Biochemistry

chejyy@ccu.edu.tw



114-2

Supplementary Material 9 : Determination of Water Hardness

115.03.22 ~ 115.03.28 (time required: 1 hour)

1. Objective

- The primary goal of this study is to understand and determine water hardness by measuring the total concentration of calcium (Ca^{2+}) and magnesium (Mg^{2+}) ions in a water sample.



2. Principle

- Water hardness is defined as the amount of calcium (Ca^{2+}) and magnesium (Mg^{2+}) ions in the water.
- It is typically expressed as parts per million (ppm) of calcium carbonate (CaCO_3), where 100 ppm represents 100 mg of CaCO_3 per liter of water (0.1 g/L).
- **100 ppm = 100 mg / L = 0.1 g / L**

EDTA

- Water hardness is determined via EDTA titration.
- EDTA (H_4Y) is slightly acidic and acts as a hexadentate ligand, using four carboxyl groups and two nitrogen lone pairs to form stable 1:1 complexes with metal cations.





Figure 1. Chemical structure of Ethylenediaminetetraacetic acid (EDTA) and its hexadentate coordination geometry with a metal cation (M^{n+}).

Indicator Mechanism of Titration

- Calmagite (H_3In) is used as the indicator. At pH 10, it forms a **burgundy** complex with Mg^{2+} .
- As EDTA is added, it displaces the indicator from the metal. At the endpoint, the free indicator turns the solution sky blue.



- Titration endpoint (at pH = 10):



3. Equipment and Chemicals

Equipment

In cabinet	Provided by TA
Beaker (100 mL) x 2	Plastic Dropper x 1
Hot plate stirrer x 1	Burette Clamp
Iron stand	Burette (50 mL) x 1
Watch glass	
Stir bar x 1	
Volumetric flask (100 mL) x 1	
Erlenmeyer flask (125 mL) x 4	
Erlenmeyer flask (250 mL) x 2	

Chemicals

6.0 M Hydrochloric acid
(HCl) ***

Ethylenediaminetetraacetic
acid (EDTA) **

0.1 % Calmagite **

Calcium Carbonate (CaCO_3) *

pH = 10.0 $\text{NH}_3/\text{NH}_4^+$ Buffer solution ***

* : Corrosive * : Toxic * : Irritant

4. Experimental Procedures

Step 1 : Preparation and Standardization of the Primary Standard (CaCO_3)

1. Initial Dissolution

- **Weighing:** Accurately weigh 0.040 g of CaCO_3 into a 100 mL beaker.
- **Solubilizing:** Add 25.0 mL of deionized (DI) water to the beaker.
- **Acidification:** Slowly add approximately 10 drops of 6.0 M HCl.
- **Reaction Control:** Cover the beaker with a watch glass and allow the solid to dissolve completely.

4. Experimental Procedures

2. Removal of Carbon Dioxide

- **Rinsing:** Use a dropper and DI water to rinse any splashes from the beaker walls.
- **Heating:** Heat the solution until it is just before boiling.
- **Degassing:** Observe the evolution of CO₂ bubbles and continue heating until the gas is expelled.
- **Cooling:** Remove from heat and allow the solution to reach room temperature.

4. Experimental Procedures

3. Quantitative Transfer and Dilution

- **Transfer:** Add 50.0 mL of DI water to the cooled solution, then carefully pour the entire contents into a 100 mL volumetric flask.
- **Bringing to Volume:** Add small increments of DI water multiple times to reach the horizontal graduation line (meniscus) of the flask
- **Homogenization:** Secure the stopper (plug) on the flask and shake or invert for at least five minutes to ensure the solution is completely uniform.

4. Experimental Procedures

4. EDTA Calibration (Titration)

- **Preparation:** Rinse the burette with 15 mL EDTA and fill it to the starting mark.
- **Aliquot:** Pipette exactly 5.0 mL of the prepared Ca^{2+} solution into a 125 mL or 250 mL Erlenmeyer flask.
- **Conditioning:** Add 6.0 mL of pH 10.0 buffer and 2 drops of Calmagite indicator.
- **Titration:** Titrate with EDTA until the solution changes from burgundy/lilac to clear sky blue.
- **Replication:** Repeat the titration twice more and average the results to calculate the final molarity of the EDTA.

- ① 0.040 g CaCO_3
- ② 25.0 mL DI water
- ③ 10 drops 6.0M $\text{HCl}_{(aq)}$




Figure 2. Explanation of Steps 1-1 ~ 1-2.




Figure 3. Explanation of Steps 1-3.

Calibrate EDTA

- ① 5.0 mL CaCO_3 (aq)
- ② 6.0 mL pH = 10.0 Buffer
- ③ 2 drops Calmagite indicator

EDTA

- Prepare three bottles of calcium ion solution
- Repeat the experiment





Figure 4. step 1-4 EDTA calibration diagram.

Step 2 : Determination of Tap Water Hardness



Using the calibrated EDTA to measure local Tap water samples.

- 1. Preparation:** Take 5.0 mL of tap water, add 6.0 mL of pH 10 buffer, and 2 drops of Calmagite indicator.
- 2. Titration:** Titrate with the standardized EDTA solution and record the volume consumed at the blue endpoint.
- 3. Replication:** Repeat the titration twice more (three times total) to calculate the average water hardness in ppm.


- ① 5.0 mL water
- ② 6.0 mL pH = 10.0 Buffer
- ③ 2 drops Calmagite indicator



Titrated with
EDTA



Titration
endpoint




- Prepare three bottles of calcium ion solution
- Repeat the experiment

Figure 5. step 2 explanation diagram.

5. Titration Processes - I

CaCO_3 is standardized using EDTA.



before titration

Figure 6. Color Change During Titration Process.



Figure 7. Titration endpoint color and titration excess color change.

5. Titration Processes - II

water hardness



Figure 8. Color change before and after titration.




Figure 9. Titration end point color.

6. Safety and Precautions (1/2)

1. **Safety:** Standard PPE (lab coat, goggles, and gloves) is mandatory at all times.
2. **Chemical Handling:** Handle all chemicals with care, as they may be toxic or hazardous. Follow all laboratory safety regulations and the instructions of the TAs at all times.
3. **Prohibited Actions:** Removal of any chemicals or laboratory equipment from the lab is strictly prohibited.
4. **Endpoint Judgment:** The EDTA endpoint blue is difficult to judge; ensure all titrations reach a similar shade to minimize error.

6. Safety and Precautions (2/2)

5. **pH Importance:** Titration must occur at pH 10 because EDTA binds Ca^{2+} more effectively at high pH, and Calmagite becomes protonated at lower pH, interfering with the color change.
6. **Solubility Risks:** If the pH is too high, Mg^{2+} reacts with OH^- to form $\text{Mg}(\text{OH})_2$ precipitate, leading to an underestimation of water hardness.

7. Experimental Data

- **Calcium carbonate is calibrated with EDTA**

1. Weight of calcium carbonate :
2. Calcium carbonate (calcium ion) mole :
3. Calcium ion concentration :
4. The number of moles of the calcium ions contained in each titration :
5. Volume of EDTA :
6. Concentration of EDTA :

- **The hardness of water is determined**

1. Volume of water :
2. Volume of EDTA :
3. The volume of EDTA required to titrate 1 liter of water :
4. The number of moles of EDTA required to titrate 1 liter of water equals the number of moles of calcium ions present in 1 liter of water :
5. The mass of calcium contained in 1 liter of water :
6. The hardness of water (1 ppm = 1 mg/L) :

The water hardness in the Minxiong Township area of Chiayi County: 207 ppm

Reference: <https://data.gov.tw/dataset/36605> (Extract day:2026/02/03)

8. Questions

1. Dissolve 0.557 g of magnesium carbonate in 6.0 M HCl(aq) and dilute to 500 mL.
 - a. Moles of MgCO_3 in the sample?
 - b. Mg^{2+} concentration?
 - c. Moles of Mg^{2+} ?

-
- 2. In this experiment, water hardness is the sum of calcium and magnesium ions.**
 - a. How can Ca^{2+} and Mg^{2+} be measured separately?**
 - b. Can this be done using the indicator from this experiment?**

**3. In Step 7, why is a pH 10.0 buffer solution added?
Why must the titration of calcium ions be carried
out under this pH condition?**

4. If the titration solution's pH is too high, Mg^{2+} reacts with OH^- to form $\text{Mg}(\text{OH})_2$ precipitate. According to solubility equilibrium, this decreases the free Mg^{2+} in solution, so less magnesium is detected during titration, leading to an underestimation of water hardness.

National Chung Cheng University
Department of Chemistry and Biochemistry

Chinese version : TA C.-H. LIAO (廖建勳)

English version: TA Y.-C. Tsai (蔡聿貞)

Both revised by: Prof. Joyce S. Yu (于淑君)

Revision date: 115.03.26