



國立中正大學

National Chung Cheng University



積極創新 修德澤人



Prof. Joyce S. Yu

Department of Chemistry and Biochemistry

chejyy@ccu.edu.tw



114-2

Expt. 4 - Quantitative Analysis of Vitamin C

115.06.05 (time required: 1.5 hours)

1. Objective

- To determine the vitamin C content in high-dose vitamin C powder and peach juice using redox titration methods.





2. Principle

- **Commercial vitamin C supplements are available in a range of dosages, such as 100 mg, 200 mg, or 500 mg. But do these products actually contain the amount of vitamin C stated on the label? This experiment aims to quantify the vitamin C content using a redox titration technique, allowing for a precise comparison between the labeled and actual values.**

- **Vitamin C, also known as ascorbic acid, is a strong reducing agent capable of reducing Fe(III) to Fe(II) or I₂ to I⁻. In this experiment, a potassium iodate (KIO₃) solution is used as the titrant in an acidic medium, where it first reacts with sodium iodide (NaI) to generate molecular iodine (I₂):**



- Iodine molecules (I_2) readily undergo a redox reaction with ascorbic acid in the solution, as illustrated in Equation 4-2:



(eq.4-2)

- Once all the ascorbic acid in the solution has reacted, any excess iodine (I_2) combines with iodide ions (I^-) to form triiodide (I_3^-). This species interacts with the pre-added starch indicator, producing a blue-black complex that signifies the titration endpoint. Using the known molar amount of iodate (IO_3^-) and its stoichiometric relationship with I_2 and ascorbic acid, the ascorbic acid content in the sample can then be accurately calculated.

$$\frac{IO_3^- (\text{mol})}{1} = \frac{I_2(\text{mol})}{3} = \frac{C_6H_8O_6(\text{mol})}{3} \quad (\text{eq.4-3})$$

3. Equipment and Chemicals

Equipment

In cabinet	Provided by TA
Graduated pipette x 1	Burette (50 mL) x 1
Safety pipette filler x 1	Burette clamp
Stir bar x 1	
Volumetric flask (100 mL) x 1	
Erlenmeyer flask (125 mL) x 2	
Beaker (100 mL) x 2	
Hot plate stirrer x 1	

Chemicals

Vitamin C Powder	Peach vitamin C juice
1 M Sodium iodide (NaI)*	0.005 M Potassium iodate (KIO ₃)**
1 M Hydrochloric acid (HCl)***	2 % starch solution*

* : Corrosive * : Toxic * : Irritative
* : Oxidizing * : Flammable

4. Experimental Procedures

A. Determination of Ascorbic Acid Content in Vitamin C Powder

1. Accurately weigh 0.1 g of vitamin C powder using an analytical balance.
2. Transfer the powder to a 100 mL beaker, add 50 mL of deionized (DI) water, and stir thoroughly until the powder is completely dissolved. Then, transfer the solution to a 100 mL volumetric flask. Rinse the beaker thoroughly with small portions of DI water and add the rinses to the flask. Finally, dilute the solution to the calibration mark with DI water and mix well.

- 3. Pipette 25.0 mL of the prepared vitamin C solution into a 125 mL Erlenmeyer flask.**
- 4. Add the following reagents to the flask:**
 - 2.0 mL of 1 M sodium iodide (NaI)**
 - 2.0 mL of 1 M hydrochloric acid (HCl)**
 - 1.0 mL of 2% starch indicator solution**
- 5. Rinse a burette thoroughly with deionized water, then rinse it three times with approximately 5 mL of 0.005 M potassium iodate (KIO_3) solution. Fill the burette with the same KIO_3 solution, remove any air bubbles, and record the initial volume (V_i) to the nearest 0.01 mL.**

6. Titrate the vitamin C solution with 0.005 M KIO_3 until a persistent dark green color appears, indicating the endpoint of the titration. Record the final volume (V_f) to the nearest 0.01 mL.




Figure 1. Solution appearance before titration (colorless)




Figure 2. Solution appearance after titration (greenish)

7. Repeat **Steps 1 through 6** to perform a duplicate titration for improved accuracy.
8. Using the titration volumes of 0.005 M KIO_3 , calculate the average ascorbic acid content in the vitamin C powder sample.

B. Determination of Ascorbic Acid Content in Juice

- 1. Prepare 150 mL of peach vitamin C juice.**
- 2. Rinse a 25 mL pipette thoroughly with deionized water, then rinse it twice with approximately 5 mL of the peach juice. Accurately pipette 50.0 mL juice into a 125 mL Erlenmeyer flask.**

- 3. Add the following reagents to the Erlenmeyer flask:**
 - 2.0 mL of 1 M sodium iodide (NaI)**
 - 2.0 mL of 1 M hydrochloric acid (HCl)**
 - 1.0 mL of 2% starch indicator solution**

- 4. Titrate the solution with 0.005 M potassium iodate (KIO₃) and record the volume used to the nearest 0.01 mL.**

5. Accurately measure another 50.0 mL of peach juice and repeat the titration to obtain a duplicate measurement.
6. Using the titration volumes of 0.005 M $\text{KIO}_3(\text{aq})$, calculate the ascorbic acid content in the peach juice, expressed in mg per 100 mL.




Figure 3. Peach juice sample before titration




Figure 4. Peach juice sample after titration (endpoint reached)

5. Precautions

- 1. Dispose of all waste liquids in the designated waste container for proper recycling.**
- 2. Shake the 2% starch indicator thoroughly before use if any precipitation is observed.**
- 3. Always wear a lab coat, goggles, and gloves during the lab session.**
- 4. Handle all chemicals with care, as they may be toxic or hazardous. Follow all laboratory safety regulations and the instructions of the TAs at all times.**
- 5. Removal of any chemicals or laboratory equipment from the lab is strictly prohibited.**

6. Experimental Data

(A) Ascorbic Acid Content in Vitamin C Powder

1. Average content of ascorbic acid in vitamin C powder: _____ (g/g)
2. Determination of ascorbic acid content in vitamin C powder:

Test	Weight of Vitamin C Powder (g)	0.005 M KIO_3 Titration Volume			Ascorbic Acid Content (g)
		V_i (mL)	V_f (mL)	ΔV (mL) $= V_f - V_i$	
1					
2					

(B) Determination of Ascorbic Acid Content in Peach Juice

1. Reported ascorbic acid content (*literature value*) : _____mg / 100 mL

2. Determination of ascorbic acid content:

Test	Volume of Peach Juice (mL)	0.005 M KIO ₃ Titration Volume			Ascorbic Acid Content (mg/100 mL)
		V _i (mL)	V _f (mL)	ΔV (mL) = V _f - V _i	
1					
2					


Average ascorbic acid content in peach juice: _____mg /100 mL

7. Experimental Error


	Experimental Value	Theoretical Value	% Error
Vitamin C Powder (g)			
Peach Juice (mg / mL)			

8. Questions

1. Vitamin C is a water-soluble vitamin. Explain this property based on its chemical structure.



2. When the four hydrogen atoms bonded to the carbon atom in methane (CH_4) are replaced with four different substituents, the resulting molecule and its mirror image become non-superimposable. Such a pair of molecules are called *enantiomers*. A carbon atom bonded to four different groups is referred to as a chiral carbon, commonly indicated as C^* . Identify the chiral carbon(s) in the structure of ascorbic acid.



3. In this experiment, Vitamin C is a good reducing agent. Explain the reason based on the chemical structure of Vitamin C.

4. Many foods contain natural reducing agents, please list five examples and identify their active reducing components.

National Chung Cheng University
Department of Chemistry and Biochemistry

Prepared : TA L.-L. Shih (李 竹)

Edited by : TA C.-H. Lee (廖晴媿)

Revised by : TA Y.-C. Tseng (曾盈臻)

and Prof. Joyce S. Yu (于淑君)

Revision date : 115.05.28